ENG

nubax

TRIO USER'S GUIDE

Please Read Thoroughly Before Use

Version 3.3

page

CONTENTS

WHAT IS THE NUBAX® TRIO?	2
THE CONTENTS OF YOUR PACKAGE	3
ASSEMBLY INSTRUCTIONS	4-10
INSTRUCTIONS FOR USE	11-13
MONITORING YOUR USE OF THE NUBAX® TRIO	14
PATIENT SPECIFIC FUNCTIONAL SCALE	15
BACK CONDITIONS AND CAUTIONS	16
FREQUENTLY ASKED QUESTIONS	17
WARRANTY	19
LIMITATION OF LIABILITY	19

WHAT IS NUBAX® TRIO?

The Nubax® Trio is a patented device that provides a gentle but strong and effective traction to the spine. Clinically proven, the Nubax® Trio provides beneficial therapy for a range of back conditions that cause pain and suffering [see Back Conditions and Cautions Pg 16]. The Nubax® Trio is also useful for maintaining a strong and healthy back as part of a normal fitness regime and can be a useful adjunct to other treatments provided by healthcare professionals.

The Nubax® Trio offers the advantage of being portable, suitable for a wide range of body sizes and shapes, lightweight and extremely easy to use. Tests conducted have shown that the traction force that the Nubax® Trio can apply to the spine is comparable to larger, more expensive machines normally found in clinics.

About this manual.

This manual accompanies the Nubax® Trio in order to provide correct assembly instructions, operational instructions and hints and tips on how to get the most from your Nubax® Trio. Please ensure that you read this manual in its entirety before commencing use of the device.

Under no circumstances should you attempt to use your Nubax® Trio if all the parts have not been supplied. If you find that there are missing parts, please contact your local distributor. Relevant contact information for users outside Australia is located on the back cover of this manual.

Nubax® recommends the regular inspection of fixtures and eyelets to ensure against wear and tear. Please also make certain that the fixtures are in tact and securely fastened before use.

If, while using your Nubax® Trio you experience severe pain and/or discomfort you should desist from using your Nubax® Trio and contact Nubax® or consult your healthcare professional for advice.

Under no circumstances should you use the Nubax® Trio in any manner other than that described by Nubax®.

THE CONTENTS OF YOUR PACKAGE


Illustration 1

Your Nubax® package should contain the following components:

Outer frame (front section)

Outer frame (mid section)

Rear leg

2 x hand grips

Inner frame with hammock

1 x chest pad with Velcro straps

Two padded shoulder pads

Waist belt

4 x locking pins

2 x long bolts with quad handles

2 x short bolts with tri-handled nuts

4 x black curved washers

1 x DVD & manual

ASSEMBLY INSTRUCTIONS

ENG

Preparation

Lay out all the components of your Nubax® Trio as shown in Illustration 1. Check the list of components and locate all of the listed components. If you have any missing or damaged components, please contact your local distributor. Relevant contact information is located on the back cover of this manual.

Step 1

Objective

Attaching Handles to inner frame

Match up the 'L' & 'R' labels on the tube ends of the hand grips with the identical 'L' & 'R' labels on the front of the inner frame. Attach hand grips to inner frame. Secure with locking pins and ensure clips are fully inserted.

page
4

ASSEMBLY INSTRUCTIONS

Step 2

Objective

Connect the 3 outer frame segments.

Match up the 'L' & 'R' labels on the two mid sections of the outer frame to the 'L' & 'R' labels on the ends of the 'U' shaped front section of the outer frame. Then insert the ends of the mid sections into the ends of the 'U' shaped section. At this stage the mid sections will move freely in the tubes. Do not be concerned about this as they will be secured in the next step.

ASSEMBLY INSTRUCTIONS


ENG

Step 3

Objective

Insert and secure the inner frame into the outer frame.


First locate the holes where the mid section and 'U' shaped section of the outer frame have been joined. Place one concave washer onto a long bolt. Push the bolt through one of the holes in the outer frame tube from the outer side, until the bolt protrudes through sufficiently for the second washer to be placed over the protruding bolt on the inside of the frame. See example picture (A) below. Repeat this step for the opposite side. When this is done, place the inner frame in the outer frame and line up the tube ends of the inner frame with the protruding bolts. Secure the inner frame into the outer frame by screwing the long bolts into the threaded ends of the inner frame tube ends.

ASSEMBLY INSTRUCTIONS

Step 4

Attach the rear leg to the outer frame and secure with retaining clips. Ensure the retaining clips are fully inserted.

ASSEMBLY INSTRUCTIONS

ENG

Step 5

Objective

To secure the chest cushion to the top of the inner frame with attached Velcro tabs.

Place chest cushion approximately mid way over the square tube on the inner frame. Fold Velcro tabs around the frame attaching the tabs to the Velcro attachments on the round tube underneath .

page
8

Step 6

Objective

Attach loops of hammock over square tube ends of inner frame.

Rotate the inner frame backwards toward the rear leg of the outer frame to a near vertical position. Then place the loops of the hammock over the ends of the square tube of the inner frame.

ASSEMBLY INSTRUCTIONS

Step 7

Objective

Attach the shoulder pads.

Match up the 'L' & 'R' labels on the ends of the shoulder pad square attachment tube with the matching 'L' & 'R' labels on the ends of the square tubes on the inner frame. Insert the shoulder pad tube end into the inner frame tube end. Depress the locating button on the under side of the square tube of the shoulder pad and push the shoulder pad into the square tube of the inner frame until the locating button locates into one of the holes in the inner tube frame.

ASSEMBLY INSTRUCTIONS

ENG

Step 8

Objective

Attach the waist belt to the main frame.

Place the waist belt across the top of the rear end of the main frame. Take a short bolt and place it through the single eyelet on one end of the belt and then placing the belt and eyelet over a hole on one side of the frame, push the bolt through the hole from the top side. Secure the bolt from underneath with a plastic handled nut. Repeat this process for the opposite side, selecting an eyelet close to the foam padding on the waist belt. (When you first go to use your Nubax Trio you should start with a short belt setting and increase the length of the belt gradually as you become familiar and confident with its use.)

Step 9

Congratulations and welcome to your Nubax Trio.

Assembly should now be complete. However, before you commence using the Nubax Trio for the first time, double check your assembly to ensure that all fastenings are tight and secure. Then read the manual and watch the DVD so that you are familiar with its correct use. Remember to check your Nubax Trio from time to time for any fastenings that may have come loose and for signs of wear and tear of the belt or fraying around the eyelets.

INSTRUCTIONS FOR USE

ADJUSTING THE NUBAX® TRIO TO FIT YOUR BODY SHAPE & SIZE

Adjusting the shoulder pads

Kneel upright into the Nubax® Trio and lift the upper section to your chest so that the shoulder pads are behind your arms and to either side of you. Note how much space there is between the side of your chest and each shoulder pad. Adjust each of the shoulder pads equally using the click buttons on the frame, closing the gap so that your chest is held firmly, but is not restricted. Be sure to adjust the shoulder pads to the same hole on each side of the frame so that your body is held symmetrically.

Adjusting the waist belt

page

II

The waist belt governs how far forward you lean, which affects the amount of traction you receive. A longer waist belt is generally required for a shorter person, while a shorter waist belt is required for a taller person. It is advisable to begin with the shortest waist belt setting and lengthen the belt one hole at a time until you have a setting that provides a stretch you feel comfortable with.* When you first use the Nubax® Trio use a short belt setting until you are familiar with the leaning forward motion. The Nubax® Trio uses your own body weight and its unique lever action to provide a gentle, yet strong gradual static traction to the spine. It is not necessary to apply additional force by pushing forward into the stretch position. Once in the device simply relax and allow your

bodyweight and the Nubax® Trio generate the stretch you require. The traction is developed by the tension between the waist belt and the shoulder pads. It is normal to experience some pressure at the back of the arms and shoulders, and across the hips and pelvis. If you are unable to find a comfortable position when you first use the Nubax® Trio, alter your position in the device by adjusting the shoulder pads and the waist belt or by adjusting where your knees are positioned until you find a position that is comfortable for you.* Follow the instructions on the following page to start using the Nubax® Trio.

How to extract yourself from a stretch

To extract yourself from the Nubax® Trio once in a stretched position:

- Place one hand on the floor, pushing up to support your bodyweight, whilst shuffling your knees forward one at a time. This will alter your centre of gravity and pull you against the waist belt, allowing you to easily raise yourself upward in the frame and into a vertical position
- Do not try to extract yourself from the frame before coming to a vertical position, as this is too awkward

As you become more familiar with the Nubax® Trio action, and more confident with its use, you will be able to simply lift yourself from a stretched position with relative ease.

* **Caution:** when you first use the Nubax® Trio you may experience some pressure in the areas around the

shoulder pads, chest pad or waist belt. This is normal. Try adjusting your position following the instructions on pages 12 and 13 until you find a position that is comfortable for you. However, if you experience severe pain or discomfort please contact Nubax® or your local distributor.

- 6 Gradually increase the duration of stretch over time. An ideal stretching time is 2-3 minutes, 2-3 times each day.

Using Nubax®

To begin using Nubax® Trio, start with the belt in the shortest position:

- 1 Place your knees well forward, right under the hammock. You should be kneeling in an upright position (figure 1)
- 2 Adjust the shoulder pads to hold you firmly but comfortably. Test how far forward you can lean. Check that your knees are well forward, underneath your hips (figure 2)
- 3 With a short belt you should not be able to lean very far forward. Return to the upright position and lengthen the belt to the next hole. Repeat this until you reach the position that applies the amount of traction you require (figure 3)
- 4 Adjusting the position of your knees forwards or backwards will also assist in finding a comfortable position that will provide you with a beneficial traction to your spine.
- 5 Gradually increase the depth of your stretch over time your stretch over time.

figure 1

figure 2

figure 3

INSTRUCTIONS FOR USE

When to use

The Nubax® Trio is a therapeutic device designed for the effective treatment for back problems and for ease of use. This product is most effective when used in the morning after moving around for at least 15 minutes, and at night immediately before going to bed.

It is also beneficial to use the Nubax® Trio immediately before playing sport or undertaking a strenuous exercise routine, as the stretch will aid the mobility of your spine. Using it again 5-10 minutes after such activity will also provide beneficial therapy.

People who work in physically demanding jobs can also benefit by using their Nubax® Trio before, after and during work.

Using your Nubax® Trio immediately prior to going to bed provides your spine with the longest period of minimal compression which aids in managing painful back conditions and maintaining a healthy flexible spine.

Healthcare professionals can use the Nubax® Trio as an adjunct to their treatments, for pre-stretching prior to physical manipulation, postural adjustments and for the patient's home treatments in between clinic visits.

Tips to achieve maximum benefit

To get the most out of your Nubax® Trio, it is necessary to relax your back muscles so that your spine can be gradually stretched. This will relieve the pressure on interdiscal tissue between your vertebrae, assist with hydration and create separation between intervertebral discs.

Take some time to experiment with the adjustments to both the waist belt and your starting position so that you can find the most comfortable position for you. Some people find it helpful to place a pillow under their knees. Always ensure that the waist belt is properly placed across your hips, at the base of the abdomen (similar to the positioning of a car seat belt).

If you feel that your shoulders are not being properly held. Adjust the shoulder pads so that they fit snugly under your armpits.

If you find discomfort from the chest pad try re-adjusting your position in the upper frame by returning to the kneeling position and moving the position of your upper torso between the shoulder pads and the chest pad against your chest. Alternatively, the chest pad can be removed and a pillow can be inserted in its place.

MONITORING YOUR USE OF THE NUBAX® TRIO

All users will find that what works for them may differ slightly from other people. It may take time to discover what works best for you. The amount of times you use your Nubax® Trio per day, the length of time you stretch and the position that you use are all unique to you. Practice until you find out what works best for you. The Nubax® Trio is also effective when used as an adjunct to therapy treatment provided by your healthcare professional.

By continuously monitoring your pain and functionality while using the Nubax® Trio, you will get an accurate indication of your improvement. Please take note of the following:

How many times a day you use the Nubax® Trio

What time of the day you use the Nubax® Trio

How many minutes you stretch for

Rate your level of pain and functionality on a scale of 0-10

Alternatively you can use the Patient Specific Functional Scale (PSFS). This is a tool to help record and monitor improvements to your functionality. Often when improving back conditions some find that although the amount of pain that they have does not decrease rapidly, their ability to partake in activities may increase and therefore quality of life is improved.

The idea of the PSFS is to write down three activities that you are unable to or are having difficulty doing as a result of your back condition. These can range from walking the dog, playing tennis or lifting a box. Rank each activity on a scale of 0 to 10 with 0 being unable to perform the activity and 10 being able to perform the activity at the same level as before the back condition or injury. Add up the scores out of 30 and keep a record. After using the Nubax® Trio for several weeks re-score the activities and see if your score total has improved. Keep reviewing your score on occasion to keep a track of any increases in functionality that you experience. See page 15 for a copy of a PSFS.

MONITORING YOUR USE OF THE NUBAX® TRIO

The Patient-Specific Functional Scale

This useful questionnaire can be used to quantify activity limitation and measure functional outcome for patients with any orthopaedic condition.

Initial Assessment:

Identify up to three important activities that you are unable to do or are having difficulty with as a result of your _____ problem.

Look at the scale below and give each activity a rating according to your ability to perform that particular activity.

Follow-up Assessments:

Please ensure that you scale each activity at regular intervals (weekly, fortnightly or monthly). In order to accurately record your progress, do not look at your previous scores prior to making a new assessment.

Patient-specific activity scoring scheme (Point to one number):

0	1	2	3	4	5	6	7	8	9	10
Unable to perform activity										Able to perform activity at the same level as before injury or problem

Activity	Initial					
1.						
2.						
3.						
4.						
5.						
Additional						
Additional						

PSFS developed by: Stratford, P., Gill, C., Westaway, M., & Binkley, J. (1995). Assessing disability and change on individual patients: a report of a patient specific measure. *Physiotherapy Canada*, 47, 258-263.

Reminder: Each time you reassess and score each activity, do so without looking at your previous scores.

Please visit the Nubax® website: www.nubax.com.au for a printable version of the PSFS

BACK CONDITIONS AND CAUTIONS

There are many causes and conditions of back pain. Using the Nubax® Trio may provide beneficial therapy to the following back conditions. Please note: The Nubax® Trio is also effective when utilised as an adjunct to other treatments and therapies as specified by your healthcare professional.

Scoliosis

When using the Nubax® Trio the shoulders and hips are correctly aligned with the spine. The Nubax® Trio allows you to have an effective means of traction in your own home, for use between clinic visits, and combined with other treatments and therapies specified by your healthcare provider amounts to an effective treatment regime for Scoliosis. The tight muscles along the concave curve of the spine are stretched and released enabling the spine to obtain a more neutral position. Once the stretch has been completed is an ideal time to then strengthen the weak longer muscles on the opposite side of the spine. These exercises can be provided by your allied health professional.

Sciatica

Research shows that the application of traction to the spine of the required force creates a distraction between the vertebra. This separation reduces pressure on the spinal nerves and back extensor muscles, which assists to reduce the pressure on the sciatic nerve, the cause of Sciatic pain. The Nubax® Trio is able to provide the necessary traction to the

spine simply by utilising the bodyweight of the user and its unique level of action.

Degeneration of the Vertebra

The application of traction to the spine of the required force creates a distraction between the vertebra. This increased space also increases the flow of fluids around the discs between the vertebra. This aids the re-hydration of compressed discs, usually associated with degenerative diseases of the spine and aids their ability to maintain themselves in a healthier state. The separation between the vertebra created by an effective amount of traction helps reduce pain and helps to maintain the separation if used regularly. This separation also allows a greater resolution of inflammation as facet joints are not constantly irritated due to their compression and that of the associated spinal nerves.

Rheumatoid Arthritis

Traction helps to increase the flow of fluids around the discs and helps to reduce pressure on the spinal column and exiting nerves. The increased fluids and reduced pressure assist in enabling re-hydration of the discs, which aids continued separation of vertebrae and further pain relief.

Prolapsed / Ruptured / Bulged Disc

Traction can assist in the reduction of pressure on injured discs, it may also be of benefit for the long term maintenance and possible reduction or alleviation of

BACK CONDITIONS AND CAUTIONS

its re-occurrence. However in the first instance of the initial occurrence of Prolapsed, Ruptured or Bulged Discs, the use of the Nubax® Trio, or any form of traction, should be considered with caution. Any inappropriate movement or excessive application of traction to an injured disc could result in additional discomfort or injury. Generally it is best to wait until the initial acute state has settled to a level that allows reasonable movement and flexibility before beginning to use traction as a method, or part of a treatment to help manage the ongoing condition. If use of the Nubax® Trio is monitored by an allied health professional it may be utilized in the acute stage of a disk injury.

page
17

Non Specific Lower Back Pain

Traction helps to reduce the pressure often felt on this section of the spine. It also helps to low to middle back muscles that often tighten throughout the course of a day as they try to support an increasingly painful back. The Nubax® Trio is effective in applying a strong, gradual, then static stretch to the lower back.

Pinched Nerve

Research shows that the application of traction to the spine of the required force creates a distraction between the vertebra. This separation of the facet joints reduces pressure on the spinal column which in turn, reduces the pressure on the exiting nerves in the back that cause pain

FREQUENTLY ASKED QUESTIONS

How often should I use my Nubax® Trio?

We recommend 2 to 3 times per day, for 2 or 3 minutes at a time. However every individual is different so experiment with frequency and length to work out what is best for you. Traction generally is most effective when used regularly and frequently, not for long periods of duration at infrequent or long intervals.

What conditions is my Nubax® Trio suitable for?

The Nubax® Trio is beneficial for :

- Osteoarthritis
- Rheumatoid Arthritis
- Ankylosing Spondylitis
- Increased Lumbar Lordosis
- Scoliosis
- Sciatica
- Disc ruptures, lesions, prolepsis
- Degeneration of the vertebra
- General Lower Back Pain
- Sacral Pain

Use of the Nubax® Trio is not recommended for users suffering from unstable Spondilolistheses, Lumbar Hypermobility or shoulder and neck pain. It is also recommended that users refrain from using the device during the latter stages pregnancy and in the early stages of post operative procedures on the back.

Is the Nubax® Trio suitable for the upper back?

The Nubax® Trio is usually beneficial from the shoulder blades (T3, T4) down to the base of the spine.

FREQUENTLY ASKED QUESTIONS

ENG

I am experiencing slight discomfort when using my Nubax® Trio. What should I do?

If the discomfort is from the external pressure being generated from either the shoulder pads or the waist belt, this is normal as the body is being held at these points and drawn apart. However, if the discomfort is more of an anatomical nature then use caution and if necessary stop using the Nubax® Trio device. It may just be a simple matter of readjusting to a more comfortable position, reducing the frequency and/or duration of use, or stop usage for a short period of time. However if the discomfort persists, or becomes increasingly painful you must cease using the Nubax® Trio device and seek assistance before proceeding with its use.

I am experiencing severe pain when using the Nubax® Trio. What should I do?

If you are experiencing severe pain you must stop using the device and seek assistance from Nubax® or your healthcare professional before proceeding with its use.

I have been using the Nubax® Trio but my existing back pain is still there. What should I do?

Sometimes a reduction in pain may not be immediate. Continue using your Nubax® Trio and monitor your pain by use of the example on page 10. Also, whilst pain reduction may not be immediate, an increase in your functionality or an increase in your ability to perform everyday tasks may occur.

Try recording this by using the Patient Functional Specific Scale on page 15.

Is the Nubax® Trio recommended and endorsed by a medical body?

The Nubax® Trio is used by Chiropractors, Physiotherapists, Osteopaths and many other healthcare professionals. It is a listed good with the Therapeutic Goods Administration of Australia

WARRANTY

Your Nubax® Trio has been designed and manufactured using the highest standards of materials and workmanship.

Nubax® guarantee to replace or repair, in whole or in part, any part or component that is proved defective through materials or workmanship for a period of up to 12 months from the date of purchase.

Warranty covers parts and labour only, and does not cover shipping costs.

The warranty is void if:

Operating instructions are not adhered to and the unit is operated incorrectly

The unit is abused and damaged through actions other than normal operations

No Proof of Purchase is provided

LIMITATION OF LIABILITY

We highly recommend that the Nubax® Trio be used in conjunction with professional advice from a qualified healthcare professional.

Your Nubax® Trio is supplied complete with an operating manual and DVD that details correct assembly and operating procedures.

The safe load limit of the Nubax® Trio is 150kg. Do not use the Nubax® Trio if your weight is greater than 150kg.

Any person using the Nubax® Trio does so at their own risk and with the understanding that Nubax® and its agents shall not be liable for any loss, injury or damages arising from use of the device.

Caution:

Any person with a pre-diagnosed condition or under current treatment should refer their use of the Nubax® Trio to their healthcare professional prior to use. The Nubax® Trio should only be used as directed.

Nubax®
PO Box 80
Leederville, WA, 6902, Australia
Ph: 1800 268 229 (Australia only)
Ph: +61 8 9314 6655 (International)
Email: info@nubax.com.au
Internet: www.nubax.com.au

United Kingdom Distributor Details:

Bioloka LTD
LAS SUITE 707, High Road
N12 0BT London
United-Kingdom

Email: info@nubax.co.uk
Ph: +443303352448

TGA Sponsor Enterprise ID Number: 30960
FDA Device Listing Number: D067518

European Regulatory Affairs Only

EMERGO EURO
Molenstraat 15
2513 BH, The Hague
The Netherlands
Ph: +31.70.345.8
Fax: 31.70.346.729

Publication Version 3.3 - 2012
S/N:PTAS/Man-B/001/0700
Copyright© 2011 Nubax®

All rights reserved, including those to produce this manual, or parts thereof, in any form without express written permission of Nubax® Pty Ltd. Pricing, availability, features, performance and specifications are subject to change without notice. Product may not exactly be as pictured in diagrams within this publication.

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Copyright © 2011 Nubax®